

# Test Case Planning

Group members	Laiba, Fariha, Brenna, Abbey, Aubrie, Jessie
---------------	--

## 1. Use Case 1 - *Log Daily Mood*

Test Case ID	TC01 - Successful Log In
Test Objective	Ensure that the system allows for a patient to login if they type in the correct credentials
Preconditions	<ul style="list-style-type: none"><li>- the app is already installed and the user is on the login screen</li><li>- the users account exists in the database</li><li>- the user has stable connection</li></ul>
Test steps	<ol style="list-style-type: none"><li>1.Enter Valid Email</li><li>2. Enter Valid Password</li><li>3. Tap Login</li></ol>
Input Values	Email: _____ password: _____
Expected results	User is successfully logged in and navigated to the mood screen

Test Case ID	TC02 - Unsuccessful Log In
Test Objective	If login fails due to an incorrect username/password, the screen should display an error message
Preconditions	<ul style="list-style-type: none"><li>-the login page is available</li><li>-the users account already exists</li></ul>
Test steps	<ol style="list-style-type: none"><li>1.Enter Valid Email</li><li>2.Enter Incorrect Password</li><li>3.Tap Login</li></ol>
Input Values	Email: _____ Password: _____ (which is incorrect)
Expected results	<ul style="list-style-type: none"><li>-Login fails</li><li>-error message displayed (email or password is incorrect)</li><li>-screen remains on the login page</li></ul>

Test Case ID	TC03 - Successful Mood Streak Update
--------------	--------------------------------------

<b>Test Objective</b>	<b>Verify that the system correctly updates the User's streak count after a successful daily mood log</b>
<b>Preconditions</b>	<ul style="list-style-type: none"> <li>- Patient has already logged in successfully</li> <li>- Patient has logged their mood for <math>n</math> consecutive days in a row.</li> <li>- Its a new day and the patient has not yet logged a mood entry for the day</li> </ul>
<b>Test steps</b>	<ol style="list-style-type: none"> <li>1. Navigate to the log mood page</li> <li>2. Select a valid mood score</li> <li>3. Tap Submit</li> <li>4. View the updated streak counter on the confirmation page or dashboard</li> </ol>
<b>Input Values</b>	Mood = 4 (valid input)
<b>Expected results</b>	<ul style="list-style-type: none"> <li>- System saves the new mood entry successfully</li> <li>- Streak counter increases by 1</li> <li>- A confirmation message appears "Mood logged successfully, you are on a 4 day streak"</li> </ul>

<b>Test Case ID</b>	<b>TC04 - Mood Streak Number Too Large</b>
<b>Test Objective</b>	<b>Validate that the mood streak counter will not increase once the threshold of four digits is reached, and the screen will display an error message.</b>
<b>Preconditions</b>	<ul style="list-style-type: none"> <li>- The user is logged in</li> <li>- The user is on the mood logger screen</li> <li>- The user's streak number is already at 999</li> </ul>
<b>Test steps</b>	<ol style="list-style-type: none"> <li>1. Select a Mood</li> <li>2. Tap Save</li> </ol>
<b>Input Values</b>	Mood = 4 (valid input)
<b>Expected results</b>	The streak number fails to update, and an error message is displayed.